


Quantification of Aerosol Species with the AMS: Background and Recent lab results relevant for OA

Presented by Manjula Canagaratna

AMS Users Meeting, Beijing, China


May 2017

How do we get species concentrations from AMS Spectra ?



- Ionization Efficiency(IE)_s = # Detected ions / Parent molecule of species s
 IE_{NO_3} , IE_{SO_4} , and IE_{NH_4} calibrated with lab aerosol: NH_4NO_3 and $(NH_4)_2SO_4$
- OA made of up 1000s of molecules, IE_{OA} can't be explicitly calibrated in lab for all
 Define IE_{OA} in terms of IE_{NO_3} and a relative OA response (Semi-empirical)

Semi-Empirical Description for IE_s




Jimenez et al., *JGR* 2003

- Slopes in this plot differ according to chemical composition (Inorganic vs. organic, variation in organic)
- Use this variation to parameterize IE in terms of primary calibrant (NO_3)

Assumptions

- IE (Y-Axis) \propto EI cross section
- # of electrons (X-Axis) $\propto MW$
- Slope $\propto IE_s/MW_s$



Eq.2)

$$\frac{IE_s}{MW_s} = RIE_s \frac{IE_{NO_3}}{MW_{NO_3}}$$

Mass Loading Calculations

$$\text{Eq. 1) } C_s = \frac{10^{12} \text{MW}_s}{\text{IE}_s Q N_A} \sum_{\text{all } i} I_{s,i}$$

$$\text{Eq. 2) } \frac{\text{IE}_s}{\text{MW}_s} = \text{RIE}_s \frac{\text{IE}_{\text{NO}_3}}{\text{MW}_{\text{NO}_3}}$$

NH₄NO₃ is primary calibrant


$$\text{Eq. 3) } C_s = \frac{10^{12} \text{MW}_{\text{NO}_3}}{\text{RIE}_s \text{IE}_{\text{NO}_3} Q N_A} \sum_{\text{all } i} I_{s,i}$$

Nitrate
Equiv.
Mass

Calibrated or Semi-Empirical

Collection Efficiency

Detection efficiency of particles that are transmitted into AMS



Huffman et al.

$$CE(d_{va}) = E_L(d_{va}) \times E_s(d_{va}) \times E_b(d_{va})$$

Lens Transmission

Overlap between
beam and vaporizer

Particle
Bounce on
oven

Mass Loading Calculations: Putting it all together

Eq. 4)
$$C_s = \frac{10^{12} \text{ MW}_{\text{NO}_3}}{\text{CE}_s \text{ RIE}_s \text{ IE}_{\text{NO}_3} Q N_A} \sum_{\text{all } i} I_{s,i}$$

- Largest Source of uncertainty for species mass concentrations

CE: Lab and Field Results

E_b is dominant contribution to CE and affected by particle bounce (depends on particle phase)

E_b affected by

References: Matthews et al., Quinn et al.

- 1) nitrate content
- 2) acidity/neutralization
- 3) Relative humidity in the sampling line
- 4) organic liquid content

Typically, ambient $CE=0.5$ (± 0.15) for all species with empirical composition dependence

(*Middlebrook et al.*)

Lab CE_s values are more variable since pure particles may not be identical in phase to complex, ambient particles (*Matthews et al., Docherty et al.*)

RIE (Relative Ionization Efficiency)

Inorganic Species

Early lab calibrations on inorganics by Frank Drewnick, Ann Middlebrook

$RIE_{NO_3} \sim 1.1$ (for ions other than m/z 30 and 46)

$RIE_{NH_4} \sim 3-5$

$RIE_{SO_4} \sim 1.2$

IE_{NO_3} , RIE_{NH_4} , are routinely calibrated for each instrument

RIE_{SO_4} can be measured for each instrument.

Complicated by non-unit CE

1) Get RIE_{NH_4} from NO_3 cal

2) Use RIE_{NH_4} from step 1 to calculate input NH_4 and SO_4 masses for sulfate calibration.

Determination of RIE for Organics

Complication: Ambient OA is made up of thousands of individual organic species. What do we use as a calibrant?

Ambient Measurements:

Usually $CE * RIE_{Org} = (0.5 * 1.4) = 0.7$

- $CE = 0.5$ assumes organic and inorganics internally mixed
- $RIE_{Org} = 1.4$ is based on unpublished laboratory work by Silva et al. (2001)

Recent Lab work : Evaluating AMS

Quantification uncertainties

- 1) What is variability in Organic RIEs across a range of O/C values?
- 2) Can simple theory be used to accurately predict measured organic RIE values?
- 3) Is RIE of individual species significantly different when in pure or mixed aerosol
(i.e. NO_3 in NH_4NO_3 vs $\text{NH}_4\text{NO}_3 / (\text{NH}_4)_2\text{SO}_4$ mix, Organics in mixed organic/nitrate particles)

IE and RIE Measurement

Uncertainties

Basic Idea: Compare measured mass with **known** input mass

Producing **known** input mass is hard, particularly for OM


- Atomizing organics (H_2O vs. other solvents)
- Calculating mass of input particles properly even with DMA size-selected particles
 - Accounting for $Q > 1$ particles in DMA output (Used single particle from pTOF)
 - asphericity (shape factors)
- loss of volatiles in lens
- Non-unit CE due to bounce

New measurements with control of input aerosol mass



Use Centrifugal Particle Mass Analyzer (CPMA) to control of input mass/particle, remove DMA $Q > 1$ ambiguity

Improved Calibration Protocol for RIE_{OA}



Data Processing

Input Mass = (CPMA Mass/particle) * CPC Number

AMS Mass = $(1/CE_s RIE_s) * \text{AMS NO}_3 \text{ Equiv. Mass}$


$CE_s * RIE_s = \text{NO}_3 \text{ Equiv. Mass} / \text{CPMA Input Mass}$

$CE_s \sim \text{AMS pTOF Single Particle Counts} / \text{CPC}$
OR
 $\text{AMS pTOF Single Particle Counts} / \text{LS}$

*Data obtained by
**Wen Xu*, Andy
Lambe* et al.***


***Manuscript in
preparation, Xu et al.***

RIE_{OA} : Lab Measurements



- Uncertainty in CE introduces uncertainty in extracted RIE s
- Some species have RIE values that are much higher than 1.4
- Use capture vaporizer to help reduce uncertainty and examine trends

RIE_{OA} vs OSc



- Capture vaporizer and Standard vaporizer show similar trends with OSc
- Lab RIE_{OA} = 4-8 for OSc < -1.5 (Ambient HOA doesn't seem to show such an obvious effect potentially since it is offset by poor lens transmission for small particles (Jimenez et al., *AS&T*, 2017))
- Lab RIE_{OA} = 1.7 for OSc > -1.5 (within 20% uncertainty of 1.4 used for ambient calcs and supports use of constant RIE value for organics in this range)

RIE_{OA} (Measurements vs. Model)

Model: $IE_S = \sigma \times J_e \times t_m$

$$RIE_S = IE_S / IE_{HNO_3}$$

σ = Ionization cross section


J_e = Electron flux

t_m = transit time of molecule through ionizer region

$\propto (1/\sqrt{MW})$


Results for OA Standards Measured in this study

- Consider t_m effect ($R^2 = 0.37$)
- Ignore t_m effect ($R^2 = 0.51$)




Poor correlation between measured and modeled RIE_{OA} values

AMS Particle Detection Process



IE_S (i.e. Ions measured/input molecules) measurements reflect

- efficiency of ionizing gas phase molecules
- thermal decomposition on vaporizer
- bounced particle vaporization off cooler surfaces
- interaction of vapors with ionizer surfaces




Included in model



processes not
accounted for in
simple model

Measured ionization efficiency reflects multiple processes
Cannot use simple model to predict RIE. Need to calibrate.

Lab Results: Mixed Inorganics



Quantification of Inorganic species not affected by matrix
 NH_4 from $(\text{NH}_4)_2\text{SO}_4$, NH_4NO_3 detected with same RIE
 SO_4 RIE constant across solutions


Lab: Binary Organic/NO₃ Particles

Can we obtain RIEs of organics by mixing with NH₄NO₃ and directly measuring org/NO₃ ratio compared to known value from DMA&CPMA size and mass selection?

Advantage:

Direct ratio to NO₃ for internally mixed particles offsets need to measure CE_B which is variable for lab particles.

Lab: Binary Organic/NO₃ Particles



Mixed Alcohol particles

- RIE_{Org} and RIE_{NH_4} agree well with measurements from pure particles

Mixed Acid particles

- RIE_{Org} and RIE_{NH_4} higher than measurements from pure particles

- Mixed Acid and nitrate particles may lose NO₃ through displacement reaction and volatilization of HNO₃
- Mixed alcohol and nitrate particles may provide an easier method of calibrating reference RIE_{Org} values

Conclusions

- Product of RIE and CE biggest source of uncertainty for calculating OA mass
- New calibration protocol allows for more accurate RIE_{OA} measurements
- RIE_{OA} vs. OSC trend examined
 - Plateau in RIE_{OA} for $Osc > -1.5$ supports use of single value. Average RIE_{OA} for this Osc range = 1.7 (within 20% uncertainty of 1.4)
 - High RIE_{OA} observed for hydrocarbons
 - Need careful evaluations of ambient OA mass closure to look for these effects
- RIE_{OA} cannot be predicted with simple models. Calibration necessary
- Mixed binary alcohol/nitrate particles may be useful method for calibrating reference RIE_{OA} values.

Acknowledgements

Wen Xu¹, Andrew Lambe¹, Phil Silva^{1†}, Weiwei Hu^{2,3}, Timothy Onasch¹, Leah Williams¹, Philip Croteau¹, Xuan Zhang^{1‡}, Lindsay Renbaum-Wolff^{1, 4}, Edward Fortner¹, Jose L. Jimenez^{2,3}, John Jayne¹, Douglas Worsnop¹, Manjula Canagaratna^{1*}

¹Center for Aerosol & Cloud Chemistry, Aerodyne Research, Inc. Billerica, MA, 01821

²Department of Chemistry and Biochemistry, University of Colorado, Boulder, CO, 80309

³Cooperative Institute for Research in Environmental Sciences, University of Colorado, Boulder, CO, 80309

⁴Department of Chemistry, Boston College, Chestnut Hill, MA, 02467

†Current address: USDA-ARS-Food Animal Environmental Systems Research Unit, Bowling Green, KY 42101

‡Current address: Atmospheric Chemistry Observations & Modeling Laboratory (ACOM), National Center for Atmospheric Research (NCAR), Boulder, CO 80307