

Overview of Experiments VI-VIII

CU- Boulder
CHEM-4181
Instrumental Analysis Laboratory

Prof. Jose-Luis Jimenez
Spring 2007

Lecture will be posted on course web page

1

Reminder of SCE Procedure and Dates

- By **10:00 AM Mon Mar 12**: email to TAs & Jose w/ groups
- By **9:00 AM Fri Mar 16**
 - email sent to TAs & Jose, explaining your idea
 - We will give you feedback about whether idea is ok, or if you need to find new one
- By **9:00 AM Wed Mar 21**
 - 2-3 page proposal sent electronically to TAs & Jose
 - 1 page introduction/motivation:
 - Question / hypothesis
 - Sample collection and storage
 - Properly cite AT LEAST 1 journal article
 - 1 page on analytical procedure & instrument
 - Properly cite AT LEAST 1 journal article
 - 1 page on chemicals needed, safety aspects, and waste generated
 - We will review the proposals and inform you of problems

2


Experiments VI-VIII

- Weeks of Mar 5 to Mar 19
- Three types instruments
- E6: PAHs in diesel exhaust
 - GC-MS
- E7: PAHs in cigarette smoke
 - HPLC w/ fluorescence & absorbance detection
- E8: Heavy metals in water
 - Electrochemistry

3

Polycyclic Aromatic Hydrocarbons

- Multiple fused benzene & pentadiene rings:



- Several are very toxic, mutagenic, and carcinogenic
- Produced in combustion processes (cars, trucks, your toaster when bread burns...) before soot

4

PAHs in Ambient Air

Model
MASPEC-13429, No. of Pages 19

ARTICLE IN PRESS



ELSEVIER

Available online at www.sciencedirect.com




International Journal of Mass Spectrometry xxx (2007) xxx–xxx



www.elsevier.com/locate/ijms

Detection of particle-phase polycyclic aromatic hydrocarbons in Mexico City using an aerosol mass spectrometer

Katja Dzepina^{a,b}, Janet Arey^c, Linsey C. Marr^{d,1}, Douglas R. Worsnop^c, Dara Salcedo^f, Qi Zhang^{a,2}, Timothy B. Onasch^c, Luisa T. Molina^{d,3}, Mario J. Molina^{d,g}, Jose L. Jimenez^{a,b,*}




E7: PAHs in Diesel Exhaust

- Summary:
 - Extract soot w/ methylene chloride
 - Solvent exchange w/ hexanes in microscale flash column
 - http://orgchem.colorado.edu/hndbksupport/colchrom/colchromprocmi_croflash.html
 - Analyze for 3 PAHs w/ GC-MS
- Practical aspects
 - Time consuming, so a person should work on standards while another starts extraction
- Safety:
 - Diesel exhaust is a “probable human carcinogen”
 - PAHs among several other toxics. Wear gloves
 - Pay attention to waste
- Quality control
 - Calibration, blank spike, and standard addition

Gas-Chromatograph (GC)


- Objective: separate compounds in mixture



Clicker Q

- What causes separation of compounds in GC? Different compounds...
 - A. have different velocities in the gas-phase (Boltzmann)
 - B. stick differently to the surface of the column
 - C. have different solubilities in helium
 - D. partition differently to the coating of the column
 - E. I don't know

Schematic of Column Chromatography



- If analyte has some affinity to the stationary phase, it will be retarded

– Equilibrium

$$K = \frac{C_s}{C_M}$$


– Kinetics

- Molecular mass transfer: diffusion

– Emerge at the detector after “retention time” t_R

9

GC Injector and Sample Loop



© 2007 Thomson Higher Education


(b)

Fig 27-4 & 27-5, Skoog

10

GC-Mass Spectrometer



- Detection by injecting sample into vacuum, creating ions, analyzing ions w/ MS
- One of several possible GC detectors (FID...)



© 2007 Thomson Higher Education

11

Concept of Mass Spectrometry



12


E6: PAHs in Cigarette Smoke

- Summary:
 - Pull cigarette smoke through a filter using vacuum
 - Filter is extracted w/ solvent
 - Analyze extract w/ HPLC w/ fluorescence & absorbance detection
 - http://www.chem.uoa.gr/Applets/AppletChrom/Applet_Chrom2.html
- Practical aspects
 - Chromatography will take 4-5 hrs. One member of group should run standards (provided) while another collects cigarette samples
- Safety:
 - Gloves: you'll be handling PAH standards, also solvents (acetonitrile, methanol, methylene chloride)
 - Remove flammable materials when lighting cigarettes
 - Pay attention to waste
- Quality control
 - 3-point calibration curves
 - Spike: inject standard and extract in the same way as sample

13

High-Performance Liquid Chrom. (HPLC)

- Conceptually very similar to GC, but several hundred atm, so more elaborate and expensive equipment



© 2007 Thomson Higher Education

Fig 28-3, Skoog 14



Clicker Q

- Why is high pressure needed in HPLC?
 - A. To remove gases from the liquid phase
 - B. To increase the solubility of the analytes in the eluent
 - C. To allow heating of the column to higher temperatures
 - D. To flow the liquid through tiny particles
 - E. I don't know

15

Absorbance Detection

- Just like UV/Vis used before, but at end of column



© 2007 Thomson Higher Education

Fig 28-7 & 28-8, Skoog 16

E8: Metal Ions in Water

- Summary:
 - Analyze H₂O samples for Pb, Cd, and Zn ions
 - Anodic Stripping Voltammetry (ASV)
- Safety
 - Concentrated nitric acid solutions
 - Goggles, bicarbonate
 - Toxic Hg and Pb
 - Gloves
- Quality control
 - Standard addition for all 3 samples

17

ASV

- Analyte is first deposited from solution into an electrode by applying a potential
- Then potential is reversed, and analyte is redissolved (“stripped”) into solution again
 - Measure potential (species) & current (amount) needed to do this




Fig 25-34, Skoog

(b) Voltammogram
© 2007 Thomson Higher Education